



Safety And Performance Information Syringe Pump

CE 0051

Table of Contents

Product Information	III
Intellectual Property	IV
Statement	IV
Maintenance Service	V
After-Sales Service Unit	V
Returns and Exchanges	V
Important Information	VI
1. Overview	1
1.1. Safety Instructions	1
1.2. Network Security	4
2. Device Description	4
2.1. Intended Use	4
2.2. Intended Use Environment	4
2.3. Intended Clinical Benefits	4
2.4. Residual Risks and Side Effects	5
2.5. Indications	5
2.6. Contraindications	5
2.7. User Qualification	5
2.8. Intended patient population	5
3. Device Installation	5
3.1. Syringe Installation	5
3.2. Method of Fixing the Syringe Pump	6
3.3. Power Supply Connection	6
3.3.1. AC Power Supply Connection	6
3.3.2. DC Power Supply Connection	6
4. Interface Function Description	6
4.1. Function Menu Interface	6
4.1.1. Record Interface	6
4.1.2. Advice System Interface	6
4.1.3. System Set Interface	7
5. Introduction to Infusion Mode	7
5.1. Speed/Time Mode	7
5.2. Trace Mode	7
5.3. Sequential Mode	7
5.4. Loading Dose Mode	8
5.5. Intermittent Mode	8
5.6. Ramp Mode	8
5.7. Dose Time Mode	8
6. Operation and Use	8
6.1. Purge	8
6.2. Start Infusion	9
6.3. Bolus	9
6.4. KOR Function	9
6.5. Shutdown	9
7. Cleaning and Disinfection	10
7.1. Daily Cleaning and Disinfection	10
7.2. Main Unit Enclosure	10
8. Maintenance and Inspection	10
8.1. Device Inspection Method	10

8.1.1. Startup Inspection	10
8.1.2. Battery Performance Inspection	11
8.2. Battery Maintenance	11
8.3. Disposal of Obsolete Device	11
9. Faults & Alarms and Troubleshooting Methods	11
9.1. Alarm Information	11
9.2. Block Alarm	12
9.3. Alarm Test	12
9.4. Alarm Mute	12
9.5. Nurse Call	12
10. Accessories	13
10.1. Recommended combination of Accessories	13
11. Product Technical Parameters	13
11.1. Anti-bolus and Blockage Alarm Delay and Possible Bolus Injection	13
11.2. Infusion Accuracy Curve	13
12. EMC	13
13. Storage and Transportation	14
14. Software Description	15
14.1.1. Portability	15

Product Information

Thank you for purchasing the syringe pump.

Please carefully review the Manual to ensure appropriate use of the device. After reviewing the Manual, please properly keep it in a convenient, easily accessible location.

This Manual provides an introduction to the products based on the most comprehensive configuration. Some information may not apply to the specific products you have purchased. If you have any questions, please contact our company.

Product Name:	Syringe pump
Model and Specification:	SP50p, SP50n
Name of Registrant:	Ambulanc (Shenzhen) Tech. Co., Ltd.
Address of Registrant:	Evergrande Fashion Huigu Building 1#101, Fulong Road, Shanghenglang Community, Dalang Street, Longhua District, Shenzhen, Guangdong 518109, China
Manufacturer:	Ambulanc (Shenzhen) Tech. Co., Ltd.
Production address:	3rd and 8th Floor, Block C, Building #5, and 1st to 10th Floor, Building #8, Skyworth Innovation Industry Park, Tangtou 1st Road, Shiyan, Baoan District, Shenzhen, Guangdong 518108, China
 EU-Representative	Shanghai International Holding Corp. GmbH(Europe) Eiffestrasse 80, 20537 Hamburg, Germany Contact pers: Qiming Cheng Telephone: +49-40-2513175 Fax: +49-40-255726
Date of manufacture	See the label on the main unit
Life time:	10 years
Software Version:	V1
Manual Version:	1.0
Revision Date of Manual:	2026.04



Note: This device is not designed for household use.

Intellectual Property


© 2024 Ambulanc (Shenzhen) Tech. Co., Ltd., All rights reserved.

The Syringe pump and its accompanying Manual are the intellectual property rights of Ambulanc (Shenzhen) Tech. Co., Ltd., including but not limited to patents, trademarks, and copyrights.

Ambulanc (Shenzhen) Tech. Co., Ltd. reserves the right to the final interpretation of this Manual.

Ambulanc (Shenzhen) Tech. Co., Ltd. reserves the right to treat the Manual as confidential information. Without the prior written approval of Ambulanc (Shenzhen) Tech. Co., Ltd., any individual or organization can not disclose or access all or part of the information in the Manual in any way.

Without the prior written approval of Ambulanc (Shenzhen) Tech. Co., Ltd., any individual or organization cannot handle the Manual in any other way, including, but not limited to, publishing, modifying, copying, distributing, leasing, recomposing, or translating.

 is the registered trademark or trademark of Ambulanc (Shenzhen) Tech. Co., Ltd. These trademarks and related Ambulanc logos are intangible property owned by Ambulanc (Shenzhen) Tech. Co., Ltd. Any other trademarks or logos referenced in the Manual are used solely for editorial purposes and remain the property of their respective owners.

Statement

Ambulanc (Shenzhen) Tech. Co., Ltd. reserves the right to modify the Manual at any time without prior notice.

Ambulanc (Shenzhen) Tech. Co., Ltd. reserves the right to modify the technology at any time without prior notice.

Ambulanc (Shenzhen) Tech. Co., Ltd. reserves the right to modify the product specifications at any time without prior notice.

Ambulanc (Shenzhen) Tech. Co., Ltd. disclaims all warranties regarding the Manual, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.

Unless otherwise specified, the terms "Ambulanc" and "the Company" in the Manual refer to Ambulanc (Shenzhen) Tech. Co., Ltd.

The product figures included in the Manual are solely for reference and may not exactly match the physical products. Refer to the physical products for the most accurate representation.

Ambulanc (Shenzhen) Tech. Co., Ltd. assumes responsibility for the product's safety, reliability, and performance only under the following conditions:

- Assembly, expansion, readjustment, improvement, and maintenance are carried out by the personnel approved by Ambulanc (Shenzhen) Tech. Co., Ltd.;
- All relevant electrical devices meet the applicable national and local requirements;
- The Syringe pump is used in accordance with the Manual.

Ambulanc (Shenzhen) Tech. Co., Ltd. assumes no responsibility for the product's safety, reliability, and performance in any of the following situations:

- Any component is disassembled, stretched, or readjusted;

- The Syringe pump is repaired or modified by personnel who are non-authorized by Ambulanc (Shenzhen) Tech. Co., Ltd.;
- The Syringe pump is not properly used according to the Manual.

Maintenance Service

Scope of Free Services:

Devices that fall within the scope of Ambulanc's warranty services.

Scope of Paid Services:

For any device that falls outside the scope of the Ambulanc warranty service regulations, our Company will provide services on a fee-based basis.

Even within the warranty period, repairs for product failures caused by the following reasons will require payment:

- Artificial damage;
- Violent or improper use;
- Grid voltage beyond the limits;
- Irresistible natural disaster;
- Use of any spare parts/consumables not approved, or any product service performed by personnel not authorized by Ambulanc (Shenzhen) Tech. Co., Ltd.



Warning

Failure to implement a satisfactory repair/maintenance plan by the individual hospital or institution using this device may result in abnormal malfunction of the device and even endanger human health.

Warranty

Manufacturing Process and Materials:

Ambulanc (Shenzhen) Tech. Co., Ltd. guarantees that, under normal use and proper maintenance, its products will be free from defects in manufacturing and materials throughout the warranty period.

After-Sales Service Unit

After-Sales Service Dept., Ambulanc (Shenzhen) Tech. Co., Ltd.

Address: 6th Floor Building #8, Skyworth Innovation Industry Park, Tangtou 1st Road, Shiyan, Baoan District, Shenzhen, Guangdong 518108, China

Service Hotline: 400-9969-120

Email: service.intl@amoulmed.com

Website: <http://www.amoulmed.com>

Returns and Exchanges

Procedure for Returns and Exchanges

In case of returning any product to Ambulanc (Shenzhen) Tech. Co., Ltd., follow the procedures below:

Obtain the right for returns and exchanges: In accordance with the Law of the People's Republic of

China on the Protection of Consumer Rights and Interests, users may contact the Customer Service Department of Ambulanc (Shenzhen) Tech. Co., Ltd. to exercise the legal rights for returns and exchanges. Please provide the Ambulanc product serial number, which can be found on the main unit enclosure or the outer packaging box. Returns and exchanges will not be accepted if the serial number is not clearly identifiable. Please specify the product model, and briefly explain the reason for the return.

Freight: The user is responsible for all freight costs associated with returning or exchanging the product to the Customer Service Department of Ambulanc (Shenzhen) Tech. Co., Ltd., including any applicable customs fees.

Important Information

1. Once the product is purchased, the customer assumes full responsibility for its maintenance and management.
2. Even if the product is still under warranty, the following situations are not covered by the warranty:
 - Damage or loss resulting from improper or rough handling.
 - Damage or loss resulting from force majeure events, such as fire, earthquake, flood, or lightning.
 - Failure to comply with the service conditions specified for the product, including any damage or loss resulting from inadequate power supply, improper installation, or unfavorable environmental conditions, may lead to system failure.
 - Damage caused by improper packaging during return shipping.
 - Damage or loss arising from the use of the product outside the territory of the original purchase.
 - Damage or loss resulting from products not purchased through Ambulanc or its authorized distributors or agents.
3. Only qualified healthcare providers with appropriate certificates can use this product.
4. Unauthorized modification of the product's software, hardware, or any other components is strictly prohibited.
5. Under no circumstances will Ambulanc be liable for any problems, damage, or loss resulting from the reinstallation, alteration, or repair of the product by personnel not authorized by Ambulanc.
6. Doctors are responsible for the clinical treatment. Ambulanc is not responsible for the clinical treatment.
7. Back up important data to external storage media, such as clinical records and notebooks.
8. Ambulanc (Shenzhen) Tech. Co., Ltd. assumes no responsibility for any loss of data stored in the product resulting from operator errors or other abnormalities.
9. Ambulanc (Shenzhen) Tech. Co., Ltd. assumes no responsibility for any damages not arising from product defects or caused by improper operation during use.
10. Ambulanc (Shenzhen) Tech. Co., Ltd. assumes no responsibility for any damage resulting from continued use of the product beyond its shelf life.
11. If the warranty claim is rejected, Ambulanc (Shenzhen) Tech. Co., Ltd. will assume no responsibility for the cost of round-trip shipping.

12. The Manual includes warnings about potential, foreseeable risks. Remain highly vigilant all the time to guard against any unspecified risks. Ambulanc (Shenzhen) Tech. Co., Ltd. assumes no responsibility for any damage or loss caused by negligence or disregard of the preventive measures specified in the Manual.

13. Once the administrator of this system changes, the Manual must be handed over.





14. A notice to the user and/or patient: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

1. Overview

1.1. Safety Instructions



Review these safety instructions carefully. These safety instructions are an integral part of the device and must always be accessible. For safety concerns, note the following:

In the Manual, the safety instructions are indicated as follows:

Prompt	Description
 Warning:	<ul style="list-style-type: none"> •Indicates potentially hazardous or unsafe operations which, if not avoided, could result in death, serious personal injury, or property damage.
 Caution:	<ul style="list-style-type: none"> •Indicates potentially hazardous or unsafe operations which, if not avoided, could result in minor personal injury, product malfunction, damage, or property loss.
 Note:	<ul style="list-style-type: none"> •Emphasizes important precautions and provides instructions or explanations for better use of this product.
 Warning:	<ul style="list-style-type: none"> • The device should be connected to a power socket with protective grounding. The use of mobile multi-position sockets is not permitted. If the power socket is not grounded, please avoid using it and power the device with batteries instead. • If the syringe pump is intended to be connected to an independent power supply, the power supply should be specified as part of the syringe pump. •When using an external power supply for the syringe pump, ensure it is connected via a plug interface, allowing for quick disconnection in the event of a failure. • Avoid using the device in an oxygen-rich environment or any area where flammable or explosive materials, such as anesthetics, are present, to prevent the risk of fire or explosion. • The device is not suitable for use in a nuclear magnetic resonance (NMR) environment. • Avoid using mobile multi-position sockets (MSOs) or AC power extension cords. Ensure that the total individual ground leakage currents do not exceed the allowable limit. •Do not open the enclosure of the device, as doing so may result in electric shock. Maintenance or upgrades to the device must be performed by maintenance personnel who have been trained and authorized by Ambulanc. • All settings should be verified and confirmed to be correct before starting the infusion. • Please arrange the power cord and other cables carefully to prevent entanglement, patient suffocation, cable interference, or electrical disruptions. • During the infusion process, blockage resulting from knotting of the infusion pipeline will lead to an increase in internal pressure within the extension pipeline. Removing the blockage at this time could lead to an excessive injection of liquid medication into the patient's body; therefore, appropriate preventive measures should be implemented. • Avoid simultaneously touching the patient and the device interfaces (USB interface and multi-functional interface) to prevent potential harm to the patient caused by leakage current. • The device is intended for use by healthcare professionals, medical electrical specialists, or professionally trained clinical medical personnel in designated settings. Personnel operating the device should be properly trained. Unauthorized or untrained personnel are strictly prohibited from performing any operations. • During the defibrillation process, do not touch the patient or their non-defibrillation devices to avoid the risk of electric shock injuries. The defibrillation process will not affect

Prompt	Description
	<p>the basic performance of the device, such as infusion, alarms, signal transmission, and other related functions.</p> <ul style="list-style-type: none"> • Regular daily cleaning and maintenance of devices is necessary. • If repairs are needed, original accessories should be used. • Utilizing different alarm settings and default configurations for the same or similar devices within the same area (such as an intensive care unit or operating room) can pose significant safety risks. • The alarm settings for different syringe pump within the same area may vary to accommodate the patient's condition. Before starting the infusion, verify that the alarm settings are appropriate for the patient receiving the infusion. • When setting the alarm sound, ensure that the volume is sufficiently loud in the current environment, allowing the operator to promptly detect the alarm when it occurs. There may be risks associated with setting the alarm sound to a lower volume. <p>The device is portable but not explosion-proof. It is not designed to be carried by patients for use during work.</p> <p>The device must be operated exclusively by medical professionals, who should be kept at an appropriate distance from the device, such as doctors or nurses, to prevent patients from operating the device on their own.</p> <p>The user should confirm that the performance provided by the pump is suitable for the intended use and is not used in any way or for any purpose other than the intended use. Failure to do so will lead to reduction of pump function, resulting in patient or user casualties.</p> <p>Any defective pump should not be used. If the syringe pump is detected to be faulty, the pump will sound an alarm and the indicator light will turn red. If this happens, please turn off the pump, disconnect the power supply and hand it over to a qualified engineer. Performance error of the syringe pump may lead to complications and patient casualties.</p> <p>The device cannot be used in situations where there are mixed gases such as flammable anesthetic gas, oxygen and ammonia oxide gas.</p> <p>Protect patients from injuries caused by overcurrent or undercurrent by setting parameters correctly and using calibrated syringes.</p> <p>Once an abnormality is found, pause should be the preferred safety measure.</p> <p>To avoid possible fault of the syringe pump, do not expose the pump to X-rays, γ-rays or ionizing radiation, or subject it to RF interference or strong electric (magnetic) field generated by diathermy apparatus or mobile phone. If the pump is to be used in or with a magnetic resonance imaging (MRI) device, it must be protected from interference by the magnetic field generated by such device. Fault of the syringe pump may lead to incorrect injection or insufficient injection, resulting in patient casualties.</p> <p>Keep away from objects that can generate strong electromagnetic waves or noise during use to avoid misoperation. Such as nuclear magnetic resonance devices, microwave generators, radiation-emitting devices (e.g., X-ray machines, CT scanners, etc.), to prevent misoperation.</p> <p>Keep the device away from high-frequency surgical instruments to avoid misoperation.</p> <p>It is prohibited to use voltages other than those specified on the device's nameplate; otherwise, it may cause damage or even fire.</p>

Prompt	Description
	<p>The battery should not be thrown into fire or exposed to heat; otherwise, it may cause the battery to leak, catch fire or even explode.</p> <p>Do not tear off the battery sheath; otherwise, it may cause explosion or hazardous chemical burns.</p> <p>Proper battery charging is essential to ensure that the pump operates on internal battery power for the designated period of time. If the battery fails to charge properly, it may lead to impaired pump functionality and potentially result in harm to patients.</p> <p>When inserting or unplugging AC power cord, be sure to hold the plug tightly and do not touch the AC plug with wet hands.</p> <p>It is recommended not to share a socket with other electrical devices.</p> <p>Do not disassemble or modify the device without proper authorization.</p> <p>The device should be checked daily. If it is not used for a long time, all functions should be confirmed to be intact before reuse.</p> <p>If any abnormality is found or if any functions are lost, please discontinue use and contact the supplier promptly. Otherwise, the manufacturer/seller assumes no responsibility for any loss, damage or injury caused.</p> <p>Avoid vibration, falling and bumping, direct sunlight or strong light.</p> <p>Avoid hot and humid air blowing directly from heaters, electric stoves or humidifiers.</p> <p>Avoid using in areas with chemical storage, excessive dust, vibration or high humidity.</p> <p>Avoid reusing or re-sterilizing disposable syringes. After use, dispose of these syringes according to appropriate instructions.</p> <p>Please work or store and transport in accordance with the working environment and the transportation and storage conditions specified in this manual.</p> <p>The characteristics of syringe pump, syringes and consumables may affect the physiological effects of drugs. Please confirm that these characteristics are consistent with the drug prescription, trumpet curve and blockage alarm time set according to the injection flow rate. Please monitor the device during injection and check the remaining liquid volume in the syringe. Do not rely solely on the device's alarm function.</p> <p>The volume in the connecting pipeline is a residual amount and will not be injected. Such extra fluid must be replenished when the syringe is first filled and the system is emptied.</p> <p>If an internal fault is detected, the device will automatically cease operation and trigger an alarm.</p> <p>The syringes and infusion lines used in conjunction with the device should have a valid medical device product registration certificate and should be properly calibrated before use. The use of inappropriate syringes and infusion lines will result in inaccurate infusion rate or dosage, compromising the therapeutic effect or causing harm to patients.</p> <p>To avoid the risk of electric shock, the power supply system of syringe pump should be properly grounded. If protective grounding is not available, please disconnect the syringe pump from the power cord, or use the device's internal battery for operation.</p> <ul style="list-style-type: none"> • When the device is connected to the same infusion port as other infusion systems, backflow may occur due to interactions between extension pipes, or the response time of the blockage alarm may be extended. Therefore, when it is necessary to connect with other infusion systems, a one-way valve should be used at the end of the extension pipeline or under the guidance of local hospital.

Prompt	Description
 Caution:	<ul style="list-style-type: none"> • Ensure a continuous power supply while using the device. In the event of a sudden power outage, certain settings may not be saved. • Given that electromagnetic fields can affect the device's performance, any device used near the device must meet EMC requirements. Mobile phones, X-ray, or MRI equipment can be sources of interference, as all emit high-intensity electromagnetic radiation. • Install or move the device properly to avoid falls, collisions, strong vibrations, or other mechanical damage. • The device should be dried immediately after exposure to rain or splashing water. • Some settings are secured with a password. Password changes should only be made by authorized personnel. If a password is required to access the relevant functions, please reach out to the appropriate personnel.
 Note:	<ul style="list-style-type: none"> • The Manual introduces the product based on the most comprehensive configurations, therefore, some of the information provided may not apply to the specific model you have purchased. • The device is equipped with a power-off storage function. After an unexpected power outage, the alarm limit settings prior to the outage can be preserved for the entire lifespan of the device. After startup, the alarm limit settings prior to power interruption can be automatically restored. • Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. • Ignoring the safety information should be considered beyond any further reasonable risk control measures.

1.2. Cybersecurity

Note:

● To ensure cybersecurity, unauthorized devices and personnel are strictly prohibited from connecting to the device to access information. Non-medical devices are also prohibited from being connected.

2. Device Description

2.1. Intended Use

The Syringe Pump is intended for use by medical staff for the purposes of controlling infusion rate and volume.

2.2. Intended Use Environment

The Syringe Pump is intended to be used in healthcare institutions and emergency transport.

2.3. Intended Clinical Benefits

The Syringe pump provides controlled and accurate infusion of medications or solutions through an intravenous routes.

2.4. Residual Risks and Side Effects

One or more risks arising from all identified hazardous conditions have been brought under control and within acceptable limits. After all risk control measures have been implemented, no new risks have been created.

The residual risks are within the acceptable range after risk control.

2.5. Indications

The Syringe Pump is for patients who need receiving various types of medications, solutions through an intravenous routes.

2.6. Contraindications

The infusion of analgesics, chemotherapy, insulin, parenteral nutrition, enteral nutrition, blood and blood products.

2.7. User Qualification

The Syringe Pump is intended to be used by trained and authorized healthcare professionals.

2.8. Intended patient population

The Syringe Pump is intended for adult, adolescents, children and infants.

3. Device Installation

3.1. Syringe Installation



Warning:

- After syringe installation, it is necessary to confirm whether it is installed in place. After successful installation, the device can automatically identify the syringe specifications. If the installation is not successful, a corresponding prompt will pop up on the screen.
- For safe operation of the syringe pump, the syringe and extension pipe must be properly installed according to the syringe installation instructions. Ensure that the syringe bezel (flange) is properly seated in the slot and the extension pipe is correctly placed on the fixing hook and infusion stand; otherwise, it may cause infusion error, and may lead to patient injury or even death.
- To avoid siphoning (free flow) of fluid in the syringe, be sure that you have clamped the extension pipe before installing or removing the syringe; ensure that the syringe is properly installed; and ensure that the syringe piston is correctly seated in the pump jaws.
- To prevent the pump from falling off the infusion stand or bed rail, ensure that the pump is securely installed, and always check the safety and stability of the installed pump. Failure to comply with this warning may result in damage to the pump and injury to the user or patient.
- Replace the syringe extension pipe according to the recommended time on the syringe extension pipe packaging or according to the hospital management requirements.
- Do not reuse the disposable accessories recommended by our company; otherwise, it may cause

performance degradation or cross infection of the syringe pump.

- Use the device within 51±5 cm above and below the patient's heart height. The smaller the height difference between the device and the patient's heart, the more accurate the pressure test in the infusion pipeline.

3.2. Method of Fixing the Syringe Pump



- Before securing the device, ensure that the stability of the fixed bracket is checked.
- When the stacking guide rails are used, up to 4 pumps can be stacked together, and when stacking, at least one clamp must be used for every two pumps.
- If multiple devices are to be stacked and mounted on the infusion stand, complete the stacking process before attaching devices to the stand.

3.3. Power Supply Connection

3.3.1. AC Power Supply Connection



- If there is any doubt about the integrity of the external protective wire during installation or wiring, battery power should be used to prevent the patient or user from possible electric shock injury.

3.3.2. DC Power Supply Connection



- Check both ends of the power connector for residual liquid medication to prevent safety hazards.
- Only power cords and batteries approved by our company can be used. Using power cords or batteries not approved by our company may cause device fault and even affect the personal safety of users.

4. Interface Function Description

4.1. Function Menu Interface

4.1.1. Record Interface



- This device can store a minimum of 30,000 historical records. When the storage capacity is reached, new records will overwrite the oldest ones.
- Power failure will not affect stored historical records.
- Cumulative volume cannot be cleared during infusion.

4.1.2. Advice System Interface



- Before entering new patient information, the existing patient data must be cleared. Upon patient clearance, all associated patient data will be permanently erased.
- Before exporting the patient information, insert the USB flash drive into the USB port on the back shell of syringe pump. If it is not connected, the syringe pump will prompt: Connect USB!

4.1.3. System Set Interface



Notes

- When the time near end, no operation time, un-infusion lock screen time, and infusion lock time are set to be off, the device will no longer execute this function, that is, there will be no related alarm.

5. Introduction to Infusion Mode

5.1. Speed/Time Mode



Notes

- The weight unit and concentration configuration can be set in the [User Maintain] interface, and the concentration parameter can be configured as either concentration or drug amt. and volume as needed.
- Method for changing the units of drug amt., dose rate and concentration in the weight mode: When selecting the units of drug amt., dose rate and concentration before infusion or during infusion pause, a corresponding unit list will pop up on the interface, and the user can select as needed.
 - ① Unit of dose rate in the weight mode: X/kg/min, X/kg/h and X/kg/24h are optional, where X represents ng, ug, mg, g, mU, U, kU, EU, mmol, mol, mcal, cal, kcal and mEq.
 - ② Units of drug amt.: ng, µg, mg, g, mU, U, kU, EU, mmol, mol, mcal, cal, kcal, and mEq.
 - ③ Units of concentration: ng/mL, µg/mL, mg/mL, g/mL, mU/mL, U/mL, kU/mL, EU/mL, mmol/mL, mol/mL, mcal/mL, cal/mL, kcal/mL, and mEq/mL.

5.2. Trace Mode



Notes

- In the trace mode, the rate must be set before starting the infusion. The vtbi and time can be selected as needed.

5.3. Sequential Mode



Notes

- In the sequential mode, when more than two sequences are set, only the rate of 1 sequence needs to be entered, and the remaining sequences can be started by entering the valid time. Only the time sequence is set as a pause sequence during the infusion process.

5.4. Loading Dose Mode



- After the device is started, it will enter the first dose phase preferably. When the device does not set the parameters of the first dose phase, the device will directly enter the main phase.

5.5. Intermittent Mode



- VTBI and maintain rate are optional settings. When the volume to be infused is not set, stop infusion after syringe emptying. When the maintain rate is not set, no infusion will be performed during the hold phase, and it will automatically switch to the infusion phase after reaching the set interval time.
- During the hold phase, users can click the [Start] button to end the hold phase and directly enter the infusion phase.

5.6. Ramp Mode



- Stab. time and stab. rate can only be obtained by calculation and cannot be entered.
- Up time and down time are optional settings. When none of them are set, the infusion pump performs infusion at a stab. rate.

5.7. Dose Time Mode



- Time can only be obtained by calculation and cannot be entered.
- Unit of dose rate in the dose time mode: X/min, X/h and X/24h are optional, where X represents ng, ug, mg, g, mU, U, kU, EU, mmol, mol, mcal, cal, kcal and mEq.

6. Operation and Use

6.1. Purge



- When the default purge rate is greater than the rate max of syringe, the purge rate is equal to the rate max of syringe.



- Before purging, please make sure that the infusion pump and extension pipe are disconnected from the patient end. If disconnected, click [Yes].
- To avoid embolism, be sure that all bubbles in the syringe and extension pipe are emptied before

any drug delivery. The syringe pump has an emptying function to assist in completing this process. The presence of air in the drugs may lead to complications and even patient casualties.

6.2. Start Infusion



Warning

- When setting the infusion parameters, be sure to check the medical advice to ensure that the infusion parameters are consistent with the medical advice.
- Ensure that the brand of syringe consumables currently used in the device is consistent with the selected brand of syringe consumables.
- Before connecting the patient, make sure that the extension pipe has been exhausted and properly connected to the device.
- When the ambient temperature is 40°C, the surface temperature of the extension pipeline may reach 43°C when it comes into contact with patient. If it gets too hot, do not contact it for a long time to avoid harming the patient.

6.3. Bolus



Warning

- When using the bolus function, make sure that the syringe pump is properly connected to the patient.

6.4. KOR Function



Notes

- When the KOR rate of syringe pump is set to 0 mL/h, the KOR function is turned off, that is, after the syringe pump completes the set volume to be infused, the KOR function will not be started. When the set KOR rate is greater than the current infusion rate of syringe pump, the syringe pump will continue to run at the current infusion flow rate.
- The KOR function runs continuously for 30 minutes. At the end of run, the syringe pump will stop running and an alarm of [KOR Completed!!!] will be triggered.
- The infusion volume of KOR will be included in the cumulative infusion volume. If the bolus function is activated during the operation of the KOR function, the device will restart the KOR function at the end of bolus.

6.5. Shutdown



Warning

- Please dispose of the used syringes and extension pipes in accordance with the relevant laws and regulations or the disposal procedures of the local hospital, and do not reuse them.



Notes

- In some special cases, if the device cannot be shut down normally, press and hold the ON/OFF button for no less than 10 seconds to force the shutdown. Forced shutdown may cause data loss or damage to the device and is not recommended.
- If the device has been connected to an external AC power supply, the AC power supply of the device will not be cut off after shutdown. To completely disconnect the power supply, please pull the power plug.

7. Cleaning and Disinfection

7.1. Daily Cleaning and Disinfection



Notes

- Only the detergents, disinfectants and treatment methods recommended by the manufacturer can be used to clean or disinfect the device or its accessories. The company does not provide any guarantee for damage or accident caused by the use of other materials or methods.
- Do not mix different disinfectants before use; otherwise, it may cause certain hazards.
- Do not pour liquid on the device or its accessories. If accidentally pouring liquid on the device or its accessories, disconnect the power supply before cleaning and drain the liquid; if the device is abnormal, please contact the professional maintenance personnel.
- After cleaning or disinfection, please check the device and its accessories. If any signs of aging or damage are found, stop using them or replace them immediately.
- Abrasive materials (such as steel wool or silver polish) or strong solvents (such as acetone or detergents containing acetone) should not be used for cleaning.
- To avoid damage to the accessories caused by improper disinfection, please disinfect the accessories only when necessary according to the local hospital's system.
- Under normal use, there is no limit to the number of times the main unit enclosure, clamp and bracket can be cleaned within the recommended life time.

7.2. Main Unit Enclosure



Notes

- Before cleaning, please disconnect the power cord of syringe pump.

8. Maintenance and Inspection

8.1. Device Inspection

8.1.1. Startup Inspection



Notes

- If the startup self-test of the device fails, a "System Error!!!" alarm will be triggered after the device is turned on. If this happens, please contact our after-sales service personnel.

8.1.2. Battery Performance Inspection



Note:

- Do not use the device for infusion during battery performance inspection, and do not interrupt charging or discharging.

8.2. Battery Maintenance



Warning:

- Do not open or remove the protective case of the battery when the device is working.
- Do not expose the battery to open flames or other heat sources. Do not dispose of it in fire.
- Do not short-circuit the battery terminals.
- Do not subject the battery to severe physical impact. Do not strike with an iron hammer.
- Please check the battery regularly for leakage and replace the leaking battery in a timely manner. In case of leakage or strange odor, please keep the battery away from the fire source to prevent the leaked electrolyte from catching fire; if the battery leaks or the leaked liquid enters your eyes, please rinse with clean water immediately and consult a doctor.
- Do not place the battery in direct sunlight or high temperature area; do not let the battery come into contact with water.
- Please store the battery in an area away from direct sunlight, high temperatures, and moisture.
- When disposing of batteries, follow local regulations.
- Do not place the battery in an unsafe environment.



Note:

- Only the battery specified by the manufacturer may be used. After replacement, please dispose of the battery according to the relevant regulations.

8.3. Disposal of Obsolete Device



Warning

- The disposal of batteries, accessories, packaging materials and device components must comply with the relevant local regulations or the hospital's waste disposal system.

9. Faults & Alarms and Troubleshooting

9.1. Alarm Information



Notes

- After the high-level alarm is triggered in the running state, the syringe pump stops running.
- After the low-level alarm is triggered, it will not affect the infusion of syringe pump, and the syringe pump continues to infuse.
- After the syringe pump emits a [Battery Depleted!!!] alarm, the pump stops infusing and the

shutdown is delayed for at least 5 minutes.

- Under the standard working conditions (ambient temperature: 20°C±2°C, rate: 5 mL/h, brightness: Gear 2, default alarm sound: Gear 6, Wi-Fi: OFF), after the first alarm of [Battery Low!], it can run for at least 30 minutes.
- When a blockage occurs, the device will stop infusion and trigger the [Block!!!] alarm.
- If there is overcurrent or undercurrent caused by partial or complete blockage of the infusion pipeline, it can be handled according to 9.6 Block Alarm. If the abnormal situation has not been eliminated, please contact the manufacturer in a timely manner.

9.2. Block Alarm



Notes

- The block pressure takes effect immediately after being set and saved, and does not change with the duration of power outage. Therefore, the blockage alarm level must be checked before injection to ensure that it meets the current injection requirements. Otherwise, the blockage alarm time may be delayed, resulting in patient casualties.
- Before the blockage is relieved, the device pumps back the drug liquid to relieve the block pressure using the method of controlling the stepping motor reversal based on the injection speed and blockage level.
- If a pipeline block alarm occurs, immediately close the rate regulator and briefly press the syringe piston release handle to reduce the pipe pressure to eliminate the possibility of the agent being injected into the patient's body. Then check whether the extension pipeline is bent or blocked or whether the syringe cap has not been removed, and eliminate the blockage before re-injection. Accidental drug injection may lead to patient casualties.

9.3. Alarm Test



Warning

- Other alarm tests required can be carried out according to the triggering causes and troubleshooting methods in "9.3 Alarm Information".
- Do not set alarm limits beyond the maximum limit. This could cause the alarm system to malfunction.

9.4. Alarm Mute



Notes

- Except for the [Battery Depleted!!!] alarm, other alarms can be paused by clicking the [Alm Mute] button.

9.5. Nurse Call



Warning

- Clinical medical staff must ensure that they are within the visible and audible range of the syringe pump, so that important alarms can be responded to quickly. If the alarms cannot be responded to quickly, it may lead to patient injury or even death.

10. Accessories

10.1. Recommended combination of Accessories



Warning

- Users should use the disposable syringe recommended by our company. If not, please contact our staff to calibrate the device again and confirm the infusion performance of the device. It can be used after confirmation; otherwise, it may affect the infusion accuracy, pressure detection and other functions of the device.
- Please be sure to use the accessories specified by our company. Using other accessories may be incompatible with the device, which may cause performance degradation or damage to the device.
- Before using the relevant accessories, please be sure to check the life time and intact packaging of the accessories. If the accessories have expired or the packaging is damaged, do not use them.

11. Product Technical Parameters

11.1. Anti-bolus and Blockage Alarm Delay and Possible Bolus Injection



Notes

- The blockage bolus volume, block alarm pressure and delay time will be affected by the test temperature, condition and pipeline length.
- When using a large-volume syringe to set a low flow rate for infusion, it may cause a greater delay in the blockage alarm. For example: When using a 50 mL syringe at the running speed of 0.1 mL/h, and setting the block pressure to the minimum and maximum, the delay in the blockage alarm may reach more than 8.5 hours and 27 hours respectively. If it is indeed necessary to infuse at a low rate, it is recommended to use a small-sized syringe to ensure the timeliness of blockage test.

11.2. Infusion Accuracy Curve



Warning

- Infusion accuracy may be affected by the environmental factors (e.g. humidity, pressure, temperature, brightness and infusion consumables used). In order to ensure the infusion accuracy of the device, the user shall set the parameters according to the parameter range specified in the manual. If beyond the range, the device cannot guarantee the infusion accuracy.

12. EMC



Notes

- SP50p and SP50n syringe pumps meet the relevant requirements of electromagnetic compatibility

in IEC 60601-1-2:2020, Chapter 11 of IEC 60601-1-12:2020, and Chapter 202 of IEC 60601-2-24:2012.

- The user should install and use the device according to the electromagnetic compatibility information provided in the accompanying documents.
- Since portable and mobile RF communication devices may affect the performance of SP50p and SP50n syringe pumps, it is necessary to keep away from such devices as smartphones and microwave ovens to avoid strong electromagnetic interference when use.
- See the accessories for details of the guideline and the manufacturer's statement.



Warning

- SP50p and SP50n syringe pumps should not be used close to or stacked on top of other devices. If it cannot be avoided, check and make sure that these syringe pumps can function normally under the configuration used.
- Class A devices are intended for use in industrial environments where EMC may be a potential challenge due to conducted and radiated disturbances from the SP50p and SP50n syringe pumps.
- Except for the cables sold as spare parts for internal components by the manufacturers of SP50p and SP50n syringe pumps, the use of accessories and cables other than those specified may result in an increase in emission or a decrease in immunity of SP50p and SP50n syringe pumps.
- This device is intended for use by healthcare professionals only. This device may lead to radio interference or disturb the operation of nearby equipment. It may be necessary to take mitigation measures, such as reorienting, repositioning or shielding the corresponding sites of SP50p and SP50n syringe pumps.
- SP50p and SP50n syringe pumps may still be interfered with by other devices even if they meet the emission requirements of the corresponding national standards. Table of emission frequency or bands, modulation types, frequency characteristics, and effective radiated power.

13. Storage and Transportation



Warning:

When the storage conditions exceed the requirements of the working environment and the device is transferred from the storage state to the use state, it should be placed in a working environment for more than 2 hours before use.

14. Software Description

14.1.1. Portability



Caution:

- The device should be disconnected from the patient before software upgrade. Please make sure that important data within the device is saved.
- Shutdown or outage should not be allowed during the upgrade of the bootstrap program, otherwise, this will lead to the breakdown of the device.



Note:

- Program upgrades can only be performed by the professional maintenance personnel.
- Before upgrading, please ensure that the version of the upgrade package is the one you need. To obtain the latest version of the upgrade package, please contact our after-sales service department.



Ambulanc(Shenzhen)Tech.Co.,Ltd.



Add: 3rd and 8th Floor, Block C, Building #5, and 1st to 10th Floor, Building #8, Skyworth Innovation Industry Park, Tangtou 1st Road, Shiyan, Baoan District, Shenzhen, Guangdong 518108, China

Tel: +86-755 26072210

Web site: www.amoulmed.com

Global Hotline: +86-400-9969-120

E-mail: info@amoulmed.com