

Global Markets



Amoul Ambulanc (Shenzhen) Tech. Co., Ltd.



Building #8, Skyworth Innovation Valley Phase II, TangTou Road No.1, ShiYan, Bao'An District, ShenZhen, China
+86 755-26072286
www.amoulmed.com
info@amoulmed.com

[Amoulmed](#) [ambulshenzhen.Tech](#) [Amoulmed](#) [Amoulmed](#) [Amoul](#)

© Copyright 2023 Ambulanc(Shenzhen) Technology Co., Ltd
SH SZAB HO-2000CE20230727-V1.10 Printed in China



HO-2000 Oxygen Concentrator

Portable oxygen concentrator
beyond your expectation

www.amoulmed.com

Amoul
We Innovate to Save Lives

Elegant design

8"

Intuitive display

High-definition 8" LED screen
Easy to adjust with touch panel

99h

Convenient Operation Timer

The smart timer can be set up to 99 hours, and it will automatically stop at the set time.


10m

Smart remote control

Allows you to set the device remotely up to 10 meters, making operation easier.



Outstanding performance



- Medical-standard oxygen supply
- Adjustable flow range 0.5-8L, personalized treatment
- Pressure Swing Adsorption technology, provide purified oxygen stably and efficiently at $\geq 90\%$ concentration
- Noise level $\leq 45\text{dB}$, sleep without disturbance
- Medicated nebulization

Accessories



1. Nebulizer

2. Mask (adult/child)

3. Nasal oxygen tube

4. Remote controller



Portable at 8kg,
easy to carry with handle



Compact design
with integrated humidifier



O₂ concentrator
and nebulizer, 2 in 1



A reliable alert system,
Escorts your safety

Alerts

- Power failure
- Low O₂ concentration
- Safety valve in humidifier
- Compressor overheating
- Screen Lock



Applications



Office workers/students



Difficult Breathing due to overweight



Patients with cardiovascular disease




Elderly persons



Patients with Respiratory disease (Asthma, COPD)



Poststroke patients



Patients with chronic disease

Specifications

Flow rate	0.5~8 L/min
O ₂ concentration	≥90%V/V(2 L/min)
Output pressure	40kPa±4kPa
Flow volume of nebulization	3L/min-10L/min
Noise level	≤45dB (A)
Dimensions	270mmx270mmx480mm
Weight	8Kg
Prefabrication time	≥30minutes