

Global Markets



Flexible Video Bronchoscope **AVL-28/39/52**

Innovative Airway Management and
Bronchoscopy Solution

Amoul® Ambulanc (Shenzhen) Tech. Co., Ltd.



📍 Building #8, Skyworth Innovation Valley Phase II, TangTou Road No.1, ShiYan, Bao'An District, ShenZhen, China
☎ +86 755-26072286
🌐 www.amoulmed.com
@ info@amoulmed.com

Amoulmed ambulshenzhen.Tech Amoulmed Amoulmed Amoul

© Copyright 2023 Ambulanc(Shenzhen) Technology Co., Ltd
S8 SZAB FVB AVLCE20231025-V1.10 Printed IN China

www.amoulmed.com

Amoul®
We Innovate to Save Lives

All-in-one design

Incorporating a touch screen for monitoring, LED as a light source, a powerful internal battery, and the ability to record data into a single device make the portable Flexible Video Bronchoscope versatile and accessible in diverse operations.





LCD High-definition Display

- 4-inch touch screen, clearer visuals
- Sensitive touch, even with gloves on
- Intelligent user interface, easy to operate



Highly Integrated, Multi-function

Multiple functions, including examination, suction, biopsy, and treatments




Easy operation, Easy observation

- A flexible bending section with a wide viewing angle ensures easy operation
- One-key press to take photos and video recording
- Rotatable monitor 270°horizontal, 120°vertically
- mini HDMI port for high-quality external display, useful for diagnosis or teaching



Easy Management

- Images and videos are recorded on TF card for easy review and transfer
- Usage frequency and time data statistic records for easy maintenance
- Powerful internal battery supports uninterrupted operation over 3 hours



Wide Range of Applications

3 different scopes are available to fit different needs in the field of airway management and bronchoscopy.

Model	AVL-3M28	AVL-3M39	AVL-3M52
Insertion tube outer diameter	2.8mm	3.9mm	5.2mm
Channel inner diameter	0mm	1.5mm	2.2mm
Length	650mm		
Angulation range	Up 180°, Down 140°		
Pixels	160,000	160,000	640,000